Nutritional Assessment Platform

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Referenties

1. Siri WE (1961). Body composition from fluid spaces and density: analysis of methods. In Brozek J, Henschel A (eds). *Techniques for measuring body composition*. Washington DC: Natl Acad Sciences/Natl Res Council pp 223-24.
2. Schutte JE, Townsend EJ, Hugg J et al. (1984). Density of lean body mass is greater in blacks than whites. *Journal of applied Physiology*, 56, 1647-1649.
3. Ortiz O, Russell M, Daley TL, Baumgartner RN, Waki M et al. (1992). Differences in skeletal muscle and bone mineral mass between black and white females and their relevance to estimates of body composition. *American Journal of Clinical Nutrition*, 55, 8-13.
4. Brozek J, Grande F, Anderson JT, Keys A (1961). Densitometric analysis of body composition: revision of some quantitative assumptions. *Ann N Y Acad Sci* 110: 113-40.
5. Lohman TG (1989). Assessment of body composition in children. *Pediatric Exercise Science*, 1, 19-30.
6. Anderson DE (2007). Reliability of air displacement plethysmography. *Journal of Strength and Conditioning research*. 21: 169-172.
7. Collins AL, McCarthy HD (2003). Evaluation of factors determining the precision of body composition measurements by air displacement plethysmography. *Eur J Clin Nutr*. 57: 770-776
8. Crapo RO, Morris AH, Clayton PD, Nixon CR (1982). Lung volumes in healthy non-smoking adults. *Bull. Eur. Physiopath. Respir*. 18, 419-425.
9. Fields DA, Hull HR, Cheline AJ, Yao M and Higgins PB (2004). Child-specific thoracic gas volume prediction equations for air-displacement plethymography. *Obes Res* 12, 1797-1804.
10. Adapted Fields-Stocks Equation (Proprietary) for children aged 2 to 11 years.
11. Stocks J, Marchard F, Kraemer R, Gutkowski P, Yishay EB, Godfrey S (1996). Plethymographic assessment of functional residual capacity and airway resistance. In stock J, Sly PD, Tepper RS, Morgen WJ (eds). *Infant respiratory function testing*. Wiley-Liss, New York, pp 191-239.
12. Fields DA, Goran MI, McCrory MA. Bodycompositionassessment via air-displacement Plethysmography in adults and children. *Am J Clin Nutr*2002;75(3):45367.
13. McCrory MA, Gomez TD, Bernauer EM, Molé PA. Evaluation of a new air displacement plethysmograph for measuring human body composition. *Med Sci Sports Exerc*. 1995 Dec;27(12):1686-91.
14. Ginde SR, Geliebter A, Rubiano F, Silva AM, Wang J, Heshka S, Heymsfield SB. Air displacement plethysmography: validation in overweight and obese subjects. *Obes Res*. 2005 Jul;13(7):1232-7.
15. Plasqui G, Soenen S, Westerterp-Plantenga, Westerteerp KR. Measurement of longitudinal changes in body composition during weight loss and maintenance in overweight and obese subjects using air-displacement plethysmography in comparison with the deuterium dilution technique. *Int J Obes* (2011) 35: 1124-1130.

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BOD POD® Body Composition Tracking System Analysis

SUBJECT INFORMATION	
NAME	[REDACTED]
AGE	57
GENDER	Female
HEIGHT	174.0 cm
ID_1	
ID_2	
ETHNICITY	General Population
OPERATOR	[REDACTED]
TEST DATE	[REDACTED]
TEST NUMBER	[REDACTED]

BODY COMPOSITION RESULT	
% FAT	35.3 %
% FAT FREE MASS	64.7 %
FAT MASS	26.456 kg
FAT FREE MASS	48.515 kg
BODY MASS	74.971 kg
BODY VOLUME	73.500 L
BODY DENSITY	1.0200 kg/L
THORACIC GAS VOLUME	3.612 L

TEST PROFILE	
DENSITY MODEL	Siri
THORACIC GAS VOLUME MODEL	Predicted

OPERATOR COMMENTS

Body Fat: A certain amount of fat is absolutely necessary for good health. Fat plays an important role in protecting internal organs, providing energy, and regulating hormones. The minimal amount of "essential fat" is approximately 3-5% for men, and 12-15% for women. If too much fat accumulates over time, health may be compromised (see table below).

Fat Free Mass: Fat free mass is everything except fat. It includes muscle, water, bone, and internal organs. Muscle is the "metabolic engine" of the body that burns calories (fat) and plays an important role in maintaining strength and energy. Healthy levels of fat-free mass contribute to physical fitness and may prevent conditions such as osteoporosis.

BOD POD Body Fat Rating Table*

**Applies to adults ages 18 and older. Based on information from the American College of Sports Medicine, the American Council on Exercise, Exercise Physiology (4th Ed.) by McArdle, Katch, and Katch, and various scientific and epidemiological studies.*

BODY FAT RATING	FEMALE	EXPLANATION
<input type="checkbox"/> Risky (high body fat)	> 40%	Ask your health care professional about how to safely modify your body composition.
<input checked="" type="checkbox"/> Excess Fat	30.1 - 40%	Indicates an excess accumulation of fat over time.
<input type="checkbox"/> Moderately Lean	22.1 - 30%	Fat level is generally acceptable for good health.
<input type="checkbox"/> Lean	18.1 - 22%	Lower body fat levels than many people. This range is generally excellent for health and longevity.
<input type="checkbox"/> Ultra Lean	15 - 18%	Fat levels often found in elite athletes.
<input type="checkbox"/> Risky (low body fat)	< 15%	Ask your health care professional about how to safely modify your body composition.

ENERGY EXPENDITURE RESULTS

Est. Resting Metabolic Rate (RMR) kcal/day	*Est. Total Energy Expenditure (TEE) kcal/day	Daily Activity Level
1359	1685	Sedentary
	2066	Low Active
	2365	Active
	2813	Very Active
	*Est. TEE = Est. RMR x Daily Activity Level	

Applies to adults ages 18 and older. Based on information from the Institute of Medicine (2002), Dietary Reference Intakes For Energy, Carbohydrate, Fiber, Fat, Fatty Acids, Cholesterol, Protein, And Amino Acids, Part I, pp93-206. Washington, D.C., National Academy of Sciences.